

وسملانا

# MANIPOL MASSAGER

LLM27



Hand grip and power  
speed regulator



Easy to operate  
and carry



4 changeable  
massage heads

Congratulations on your Lifelong Purchase. Our name stands for high-quality and thoroughly tested products. Please read these instructions carefully before use.

## HOW TO USE

The Lifelong Manipol Massager can give you an effective full body massage in the comfort of your own home. The sturdy design and convenient handle helps you massage with ease.

Step 1: Put the plug into the socket.

Step 2: Choose between the massage head options (rolling, scraping, wave) depending on your personal preference

Step 3: Turn the round knob on the device to turn it on as set massage speed

Step 4: Guide the device to the area of the body to be massaged. You can stand, sit or lie down and use the massager

Step 5: Glide the speed knob and adjust the speed as per your requirement

Step 6: Use on a particular area only for 15 minutes and then switch to another.

Step 7: Unplug the device and let it cool down use. Keep the device dry.



Arm



Shoulder



Back



Leg



Foot

## FEATURES

1. The comfortable and convenient handle, massages various parts of the body.
2. The massage strength is adjustable.
3. The vibration of the massager can help in breaking down fat deposits on the body, thus helping weightloss.
4. The machine has three changeable massage heads, for different results.
5. Light weight and convenient to use.

**Warning:-** Please ensure the product is covered with a protective cloth when in use to avoid injury

## PRODUCT PART INTRODUCTION



## MASSAGE HEADS



Wave  
massage head



Rolling  
massage head



Scraping  
massage head



Protective  
clothing

## Consumer Warranty Card

Dear Customer,

Thank you for choosing a Lifelong consumer product. All Lifelong Consumer products are designed and manufactured to the highest standards to deliver high quality performance, as well as easy installation and use. At Lifelong, we believe in providing not only service, but adding value to your purchase. The warranty has therefore been designed especially for you with your interest at heart.

### Warranty Service

All Lifelong consumer products are covered against manufacturing defects from their date of purchase.

Name of the product: \_\_\_\_\_

Model: \_\_\_\_\_

Warranty coverage: 6 months

Please note:

Purchase receipt is necessary for warranty verification.

Customer Care: [customercare@lifelongindia.com](mailto:customercare@lifelongindia.com)

### Customer Details

Name: \_\_\_\_\_

Address: \_\_\_\_\_

Home Number: \_\_\_\_\_

Office Number: \_\_\_\_\_

E-mail Address: \_\_\_\_\_

### Product Details

Model No.: \_\_\_\_\_

Serial No.: \_\_\_\_\_

Purchase Date: \_\_\_\_\_

Invoice Number: \_\_\_\_\_

Online Site: \_\_\_\_\_

Please log onto to [www.lifelongindiaonline.com](http://www.lifelongindiaonline.com) and complete the online warranty form with your personal and product details within 14 days.

## IMPORTANT

1. After using and before cleaning, you must unplug
2. If the machine operation is informal, please cut off power immediately, as soon as possible to contact the dealer
3. If the supply cord is damaged, please do not use
4. Pregnant women or people with skin disease, heart disease and other diseases or inflammation of the skin, please refrain from using
5. Children or the disabled should not use the machine
6. Do not use it over 15 minutes in the same position
7. Do not use this machine in damp places, such as bath room
8. It is suggested that the machine is at least an hour after means
9. If replacing the massage head, please make sure the machine has stopped.
10. After using, place the product in a clean place with no direct sunlight.
11. Do not pull the cord or pull the product with the cord.
12. Do not put water into the product.
13. Please use a soft cloth to clean the machine.



Use it Protective clothing correctly to protect the hair.



If use it wrongly, the hair will be tangled.

## TECHNICAL PARAMETERS

Voltage: AC 220V - 240V

Frequency: 50Hz

Rated power: 25W

### Terms and Conditions

1. This warranty is void if:
  - a. The completely filled warranty card is not presented at the time of servicing the product.
  - b. The product is not operated according to the instructions given under the user manual.
  - c. Damages are caused by lightning, abnormal voltage, water or other liquid intrusion, fire, flood, accident, negligence or improper handling.
  - d. Product has been damaged due to installation, repairs, alterations or modifications by unauthorized service organizations or persons.
  - e. Product label specifying the model number, serial number and production code has been removed and altered.
  - f. Defects or parts requiring replacement due to ordinary wear and tear, corrosion, rust or stains, scratches, dents on the casing or paintwork of the product.
  - g. Claims damaged and/or missing parts (accessories) after 7 days from the original date of goods received.
  - h. Defects or faults in product which have been used for commercial/industrial purposes or which have been rented/leased or have been otherwise subject to non-household/non-domestic use.
2. Repairs or replacements will be carried out by authorized service provider.
3. During the limited warranty period, Lifelong or its authorized service provider will repair without charge the defective unit inclusive of labor and parts and restore the unit to its optimum working condition. All defective parts used for the warranty repair should be surrendered to Lifelong and/or its authorised servicer.
4. All expenses incurred in collecting the unit (s) or part(s) thereof from authorised service provider as well as any other expenses and incidentals will be borne by the consumer.
5. Lifelong obligations are limited to the repair and replacement of the defective product. Except as set forth above, there are no other express or implied warranty and all warranted, conditions or other terms implied by statute or common law (including any warranty of satisfactory quality, merchantability or fitness for a particular purpose) are excluded to the fullest extent permitted by the law.
6. Lifelong total liability for damages relating to or arising out of the purchase or use of the product regardless of the type or cause of such damage of the form of characterization of the claim asserted (e.g. contract or tort) shall not exceed the original purchase price paid for the product.
7. However in no event shall Lifelong and Lifelong authorized distributors be liable for any punitive, special incidental, indirect or consequential losses or damages whatsoever (including without limitation, damages for the loss of revenue, business, profits, goodwill, or contracts, business interruptions, loss of business information or any other pecuniary loss.) Whether or not Lifelong has been advised of the possibility of such damages. These limitations shall apply notwithstanding the failure of the essential purpose of any limited warranty. This limited warranty does not affect consumer's statutory rights under the law.
8. No carrier, dealer or employee is thereof authorized to make modifications to this warranty and you should not rely on any such representation. Lifelong reserves the right to amend the terms and conditions if necessary.

Manufactured for & Marketed by :  
**Lifelong Online Retail Private Limited**  
Made in China

For queries and complaints please contact:  
customer@lifelongindia.com  
www.lifelongindiaonline.com