

وسملانا

RECHARGABLE BODY MASSAGER

LLM45



Rechargeable



3 Interchangeable
massage heads



Non slip handle

Congratulations on your Lifelong Purchase! Please go through the instructions carefully to help you use this product in a safe and effective manner.

HOW TO USE

Step 1: Ensure that the device is switched off before you plug in.

Step 2: Choose between 4 different massage heads

Step 3: Switch on the device pressing the power button. The LED indicator will light up.

Step 4: Press the mode button to change the rhythm of the massage

Step 5: Press the speed button to increase/decrease speed

Step 6: Guide the device to the area of the body to be massaged. You can stand, sit or lie down and use the massager.

Step 7: Use on a particular area only for 15 minutes and then switch to another.

Step 8: Unplug the device and let it cool down use. Keep the device dry.

SPECIFICATIONS:

Voltage: 8.4V

Power: 16W

Input polarity: 



Arm



Shoulder



Back



Thigh



Foot

Consumer Warranty Card

Dear Customer,

Thank you for choosing a Lifelong consumer product. All Lifelong Consumer products are designed and manufactured to the highest standards to deliver high quality performance, as well as easy installation and use. At Lifelong, we believe in providing not only service, but adding value to your purchase. The warranty has therefore been designed especially for you with your interest at heart.

Warranty Service

All Lifelong consumer products are covered against manufacturing defects from them date of purchase.

Name of the product: _____

Model: _____

Warranty coverage: 1 Year

Please note: Purchase receipt is necessary for warranty verification.

Customer Care: customercare@lifelongindia.com

Customer Details

Name: _____

Address: _____

Home Number: _____

Office Number: _____

E-mail Address: _____

Product Details

Model No.: _____

Serial No.: _____

Purchase Date: _____

Invoice Number: _____

Online Site: _____

Please log onto to www.lifelongindiaonline.com and complete the online warranty form with your personal and product details within 14 days.

USE INSTRUCTIONS

Charging method

The product should be fully charged if using for the first time or after a long period.

1. Insert the charging plug into the charging jack, making sure it's off.
2. The product takes approximately than 3.5 hours to charge. (When the product is fully charged, the indicator changes from red to green).
3. After the charging is completed, please pull out the charger and the charging plug.

A fully charged massager can be used for about 80 minutes.

This product cannot be used when charging, please disconnect the charger before use.

FUNCTION AND FEATURES

1. Rechargeable, easy to use and carry
2. Microcomputer control
3. 4 kinds of massage heads and a combination of various massage rhythms and speeds
4. Built in indicator light

Terms and Conditions:

1. This warranty is void if:
 - a. The completely filled warranty card is not presented at the time of servicing the product.
 - b. The product is not operated according to the instructions given under the user manual.
 - c. Damages are caused by lightning, abnormal voltage, water or other liquid intrusion, fire, flood, accident, negligence or improper handling.
 - d. Product has been damaged due to installation, repairs, alterations or modifications by unauthorized service organizations or persons.
 - e. Product label specifying the model number, serial number and production code has been removed and altered.
 - f. Defects or parts requiring replacement due to ordinary wear and tear, corrosion, rust or stains, scratches, dents on the casing or paintwork of the product.
 - g. Claims damaged and/or missing parts (accessories) after 7 days from the original date of goods received.
 - h. Defects or faults in product which have been used for commercial/industrial purposes or which have been rented/leased or have been otherwise subject to non-household/non-domestic use.
2. Repairs or replacements will be carried out by authorized service provider.
3. During the limited warranty period, Lifelong or its authorized service provider will repair without charge the defective unit inclusive of labor and parts and restore the unit to its optimum working condition. All defective parts used for the warranty repair should be surrendered to Lifelong and/or its authorised servicer.
4. All expenses incurred in collecting the unit (s) or part(s) thereof from authorised service provider as well as any other expenses and incidentals will be borne by the consumer.
5. Lifelong obligations are limited to the repair and replacement of the defective product. Except as set forth above, there are no other express or implied warranty and all warranted, conditions or other terms implied by statute or common law (including any warranty of satisfactory quality, merchantability or fitness for a particular purpose) are excluded to the fullest extent permitted by the law.
6. Lifelong total liability for damages relating to or arising out of the purchase or use of the product regardless of the type or cause of such damage of the form of characterization of the claim asserted (e.g. contract or tort) shall not exceed the original purchase price paid for the product.
7. However in no event shall Lifelong and Lifelong authorized distributors be liable for any punitive, special incidental, indirect or consequential losses or damages whatsoever (including without limitation, damages for the loss of revenue, business, profits, goodwill, or contracts, business interruptions, loss of business information or any other pecuniary loss.) Whether or not Lifelong has been advised of the possibility of such damages. These limitations shall apply not withstanding the failure of the essential purpose of any limited warranty. This limited warranty does not affect consumer's statutory rights under the law.
8. No carrier, dealer or employee is thereof authorized to make modifications to this warranty and you should not rely on any such representation. Lifelong reserves the right to amend the terms and conditions if necessary.

Manufactured for & Marketed by :
Lifelong Online Retail Private Limited

Made in China

For queries and complaints:
please contact: customercare@lifelongindia.com
www.lifelongindiaonline.com