

وسملینا

# FOOT MASSAGER

LLM90



Kneading, rolling and  
vibrating actions



Preset and  
manual modes



Washable  
foot cloth

Dear Customer,

Congratulations on your Lifelong Purchase. Our name stands for high-quality and thoroughly tested products. Please read these instructions carefully before use.

## HOW TO USE

STEP 1: Ensure that the device is switched off before plug in

STEP 2: Put your feet in the black cloth area

STEP 3: Switch on the device by sliding the switch to ON position

STEP 4: This massager has 4 automatic preset levels P1, P2, P3 & P4.

STEP 5: Click on the manual button and you can change the speed and directions as per your requirement

STEP 6: Press the Power Button Switch off the machine



## CLEANING

To clean the unplugged device after it has cooled down, wipe with a slightly dry cloth. However, ensure that water never gets inside the device. The foot cloth can be removed from the zipper and washed separately

## IMPORTANT

1. Ensure that the device is switched off before you plug in.
2. Do not use the object with wet hands to plug in the power button
3. Do not use the device if it is in a poor condition
4. Do not use the device with wet feet
5. Never operate the device unsupervised, especially if there are children in the vicinity.
6. Do not use with animals or pets.
7. Persons with cardiac pacemakers are requested to consult their doctor before using the massager device
8. You must never fall asleep during use.
9. Connect the device only to the voltage stated on the type nameplate.
10. As with every electrical device, this massage device should be used carefully and cautiously in order to avoid the risk of electric shock:

## Consumer Warranty Card

Dear Customer,

Thank you for choosing a Lifelong consumer product. All Lifelong Consumer products are designed and manufactured to the highest standards to deliver high quality performance, as well as easy installation and use. At Lifelong, we believe in providing not only service, but adding value to your purchase. The warranty has therefore been designed especially for you with your interest at heart.

### Warranty Service

All Lifelong consumer products are covered against manufacturing defects from the date of purchase.

Name of the product: \_\_\_\_\_

Model: \_\_\_\_\_

Warranty coverage: 1 year

Please note:

Purchase receipt is necessary for warranty verification.

Customer Care: [customercare@lifelongindia.com](mailto:customercare@lifelongindia.com)

### Customer Details

Name: \_\_\_\_\_

Address: \_\_\_\_\_

Home Number: \_\_\_\_\_

Office Number: \_\_\_\_\_

E-mail Address: \_\_\_\_\_

### Product Details

Model No.: \_\_\_\_\_

Serial No.: \_\_\_\_\_

Purchase Date: \_\_\_\_\_

Invoice Number: \_\_\_\_\_

Online Site: \_\_\_\_\_

Please log onto to [www.lifelongindiaonline.com](http://www.lifelongindiaonline.com) and complete the online warranty form with your personal and product details within 14 days.

11. Unplug the mains cable immediately after use.

12. Prevent any contact with water (except when cleaning with a slightly damp cloth). Water must never get inside the device. Never submerge the device in water. Never use the device in the bath, in the shower, in a swimming pool or over a basin filled with water. However, if water has got inside the casing, unplug the device immediately and contact customer service.

13. The device and in particular the mains cable should be checked for any damage before each use. Never plug the device into the mains power supply if the cable is damaged.

14. Never hold or carry the device by mains cable.

15. If a fault should occur during use, switch off the device immediately and unplug it from the socket.

16. The device is intended only for the purpose described in these instructions

for use. The manufacturer can accept no liability for injury arising from incorrect or careless use.

## Terms and Conditions

1. This warranty is void if:
  - a. The completely filled warranty card is not presented at the time of servicing the product.
  - b. The product is not operated according to the instructions given under the user manual.
  - c. Damages are caused by lightning, abnormal voltage, water or other liquid intrusion, fire, flood, accident, negligence or improper handling.
  - d. Product has been damaged due to installation, repairs, alterations or modifications by unauthorized service organizations or persons.
  - e. Product label specifying the model number, serial number and production code has been removed and altered.
  - f. Defects or parts requiring replacement due to ordinary wear and tear, corrosion, rust or stains, scratches, dents on the casing or paintwork of the product.
  - g. Claims damaged and/or missing parts (accessories) after 7 days from the original date of goods received.
  - h. Defects or faults in product which have been used for commercial/industrial purposes or which have been rented/leased or have been otherwise subject to non-household/non-domestic use.
2. Repairs or replacements will be carried out by authorized service provider.
3. During the limited warranty period, Lifelong or its authorized service provider will repair without charge the defective unit inclusive of labor and parts and restore the unit to its optimum working condition. All defective parts used for the warranty repair should be surrendered to Lifelong and/or its authorised servicer.
4. All expenses incurred in collecting the unit (s) or part(s) thereof from authorised service provider as well as any other expenses and incidentals will be borne by the consumer.
5. Lifelong obligations are limited to the repair and replacement of the defective product. Except as set forth above, there are no other express or implied warranty and all warranted, conditions or other terms implied by statute or common law (including any warranty of satisfactory quality, merchantability or fitness for a particular purpose) are excluded to the fullest extent permitted by the law.
6. Lifelong total liability for damages relating to or arising out of the purchase or use of the product regardless of the type or cause of such damage of the form of characterization of the claim asserted (e.g. contract or tort) shall not exceed the original purchase price paid for the product.
7. However in no event shall Lifelong and Lifelong authorized distributors be liable for any punitive, special incidental, indirect or consequential losses or damages whatsoever (including without limitation, damages for the loss of revenue, business, profits, goodwill, or contracts, business interruptions, loss of business information or any other pecuniary loss.) Whether or not Lifelong has been advised of the possibility of such damages. These limitations shall apply notwithstanding the failure of the essential purpose of any limited warranty. This limited warranty does not affect consumer's statutory rights under the law.
8. No carrier, dealer or employee is thereof authorized to make modifications to this warranty and you should not rely on any such representation. Lifelong reserves the right to amend the terms and conditions if necessary.

Manufactured for & Marketed by :  
**Lifelong Online Retail Private Limited**

Made in China

For queries and complaints please contact:

customercare@lifelongindia.com

www.lifelongindiaonline.com