

Manufactured Market by:  
**Lifelong Online Retail Private Limited**

For queries and complaints:  
please contact: [customercare@lifelongindia.com](mailto:customercare@lifelongindia.com)  
[www.lifelongindiaonline.com](http://www.lifelongindiaonline.com)

lifelong

## Electronic Pulse Tens Massager



Accelerate  
the blood  
circulation



Relax the  
regional  
muscle



Eliminate  
pain



Eliminate  
swelling



Eliminate  
inflammation



Eliminate  
fatigue

**1 YEAR  
WARRANTY**

## INTRODUCTION:

Welcome to the Lifelong family,

Dear Customer,

Congratulations on purchasing Lifelong Tens Massager and a warm welcome to the Lifelong family. At Lifelong, we take pride in developing quality products for your home, health and happiness. Through constant innovation in our products, we have been able to develop a wide spectrum of products across Home, Kitchen, Grooming and Wellness categories. We strive to enrich the lives of the modern consumer every day. Keeping that in mind, we have offered all our products at a fair and honest price.

We hope that the Tens Massager serves you and your family's needs without any hassles. Thanks again for becoming a part of the Lifelong family. Enjoy using the product!

Warm Regards,  
Team Lifelong

## PRODUCT PREVIEW:

Electronic Pulse Massager is a new generation of health care of high-tech products based on Meridian, bio-electricity, physics, combined with long-term experience in clinical practice.

It uses the microcomputer control technology to program the weak current regularly that is acting on the body for stimulating, and realistically simulates acupuncture, cupping, massage, thump, effective activating human cells, making tissue cells excited to promote blood circulation, thus dispels pains and prolongs life.

## TREATMENT MECHANISM:

Modern medicine has proved, Meridian is a sensing and low-resistance conductor. In traditional acupuncture. Every twisting or pull the stick once, which is to generate a stimulating electrical pulse. In addition, in human body life, some particular current is happened all the time. These weak current called "bio-electricity" are sent from the heart, brain, muscles, nerves and play an important role for the normal operation of body. Throughout the process of human body life, it can't be separated from the the bioelectronic movement, it always maintains a relative balance in imbalance. Viewed from the meridian point, the imbalance of the food chain will lead to the meridian disorder, causing various human discomforts or pains. Electronic Pulse Massager makes the weak low-frequency pulse current to be applied to the human body, using it to stimulate the body to produce reaction, thus expand blood vessels, promote blood circulation, improve tissue viability, accelerated metabolic waste and efflux of inflammatory substances, reaching to anti-inflammatory, consumer swollen. The high-energy current can cause the masking effect, though the interference and the shutter mechanism to generate the effect of analgesia timely.

While the fluid circulation can stimulate an endorphin content increased in the morphine-like substances in serum, producing sustained analgesic effects; high energy can stimulates nerve excitability and muscle contraction caused by neuromuscular to produce movement effects and can eliminate fatigue and treatment of peripheral nerve damage and paralysis.

## FEATURES:

1. Porcelain appearance, easy to carry
2. Micro-computerized system, LCD display, easy to operate.
3. One automatic compound massage mode with 4 different manual massage modes. Applicable to different needs of people.
4. Wonderful combination of electrical pulses, allowing you to truly appreciate the acupuncture, scraping, cupping, massaging, hammering.
5. Battery power display.

## PRODUCT APPLICABLE SCALE:

1. Limb pain disorders: cervical spondylosis, bone hyperplasia, lumbar muscle strain, lumbar disc herniation, sciatica, femoral head necrosis, rheumatism and rheumatoid arthritis, frozen shoulder, sprained, old injuries, leg pain.
2. Nervous system diseases: stroke buckle sequelae, hemiplegia, numbness in hands and feet, trigeminal neuralgia, facial paralysis, varicose veins, hypertension, neurasthenia, headache, insomnia.

3. Class of diseases of the gastrointestinal tract: stomach problems, gastritis, gastrointestinal ulcers, gastrointestinal spasms, stomach and abdominal pain, constipation.
4. Cosmetic effect: weight loss breast clear acne.
5. Other: physical discomfort, low immunity, general fatigue, toothache, menstrual disorders, impotence, sexual function decline.

## PRODUCT APPLICABLE SCALE:

1. Putting two of AAA battery into control panel.
2. Tearing the protective film of electronic conduction patches, then stick it onto sore spot. Press the "0" button on the control panel to operate.
3. Program of work start, enter into automatic mode, namely below: A:Automatic massage mode. B:First massage mode. C:Time on 15 minutes.If you want to stop massage, please press the "D" button again.
4. The usage of "g 9" button "g 9" button for choosing five modes, namely automatic, cupping, kneading, pressing, massage, the automatic mode is cycle cupping, kneading, pressing, massaging.
5. Instensity Operate Start that is set to 1 intensity. Instensity has sixteen degrees. Press the "+" key to strengthen intensity, and press the "-" key to weaken intensity. Slowly adjust the massage intensity and it is advisable to self accept for intensity, the left and right side can be used simultaneously, and also may be used alone side.

6. Operate for timer Auto stopped after 15 minutes of operation. Please press the "D" button for use once more. Recommended 15 minutes every day.



## BATTERY LIFE AND IMPORTANT ITEMS:

1. "AAA LR03" batteries are used about 3 months.(1 times per day, every 15 minutes).
2. When LCD display shows " Q ",which means the battery is hord of supply,please replace the new battery.
3. Take off battery for long time no use to avoid leakage.
4. Do not use the battery without designation.
5. Do not insert the batteries with their polarities in the wrong direction,please read the symbol " +" and"—" clearly.

6. Please cuff off the power before replace the batteries.
7. Battery waste is processed in accordance with city regulations on environmental protection.

## METHOD OF CONDUCTIVE FILM:

1. Conductive film surface should be kept clean to avoid glue dust, oily substance, goo dirt. Otherwise viscous will decline, when the conductive film surface is dirty, cleared with water, and reuse when it dry.
2. Conductive film is generally used about 80 times (every 15 minutes).When multiple cleaning conductive film does not restore its viscous, please consider to purchase a new conductive viscose electrode to the dealer or manufacturer.

## CLEARING MAINTENANCE AND CUSTODY:

### 1. Clearing and Maintenance




1. If the machine is dirty, using the soft cloth or dry cloth to wipe.
2. Conductive film surface should be kept clean to avoid glue dust, oily substance, goo dirt. Otherwise viscous will decline, when the conductive film surface is dirty, cleared with water, and reuse when it dry.




## 2. Custody



1. Please place it where without children' reach.
2. Please do not put in the place of direct sunlight, hot temperature & humid place.
3. Please keep it in dry&ventilate place.
4. Please do not dispart the instrument, and not repair of modification, or causing the breakdown.

## NOTES ON SAFETY:

The icon in the introduction is to make you safe and right using it. Icon,picture and meaning as following own:

Alarm mark	Mean
 Dangerous	It is danger to have accident with incorrect operation
 Alarm	It is danger to have accident with incorrect operation
 Notice	It is possible to have accident or damage with incorrect operation.

Picture mark	
	is forbidden to do anything, left icon is prohibition of disassembly.
	is obliged to comply.
	type BF applied part.

 <b>Dangerous</b>	
<ol style="list-style-type: none"><li>1. Not use it with other electrical medical instrument at the same time, such as the pace maker etc,or would cause danger.</li><li>2. If the patient use the high frequency operate instrument with our therapy instruemnt,then would cause burning of the massage place, also broke the instruemnt.The output of therapy instrument would unstable if having low wave or micro wave instrument near in 1 meter when using it.</li><li>3. Closing to the chest electrode will increase the risk of cardiac fibrillation.</li></ol>	

⚠ Alarm	
Pregnant, sensitive patient, Cardiovascular patients or other patient who accept the treatment of doctor should consult the doctor before using it.	!
1. Do not use beside heart, head, mouth or private part, skin disease. 2. Do not use in the high temperature, flammability and radiant place.	⊘
3. Do not use when bathing and sleeping 4. Do not use the therapy instrument during pregnant or Menstrual period	
Do not disassemble the instrument, and not repair or modification, or causing the breakdown.	⚡⊘

⚠ Notice	
1. Do not move the instrument when in use. 2. Children should use with adults' guide. 3. Any abnormal discomfort because of it, stop using and ask doctor. 4. Taking off battery for long time no use to avoid leakage.	

## TECHNICAL SPECIFICATION:

1. Operating condition
<ul style="list-style-type: none"> <li>• Temperature: +5°C~+40°C</li> <li>• Relative humidity: ≤80%</li> <li>• Atmospheric pressure: 86kPa~106kPa</li> <li>• Adaptor: DC 3V (2x AA LR03 battery) 20 mA</li> </ul>
2. Electric performance :
<ul style="list-style-type: none"> <li>• Pulse frequency: 1Hz~1000Hz</li> <li>• Pulse extent: 20μs~400μs</li> <li>• Impedance: 500Ω, accuracy: ±10%</li> <li>• Timing: 15 minutes, accuracy: ±10%</li> </ul>
3. Transportation & Custody :
<ul style="list-style-type: none"> <li>• Circumstance temperature: -20°C~+55°C</li> <li>• Relative humidity: ≤93%</li> <li>• Atmospheric pressure: 50kPa~106kPa</li> </ul>

## Consumer Warranty Card

Dear Customer,

Thank you for choosing a Lifelong consumer product. All Lifelong Consumer products are designed and manufactured to the highest standards to deliver high quality performance, as well as easy installation and use. At Lifelong, we believe in providing not only service, but adding value to your purchase. The warranty has therefore been designed especially for you with your interest at heart.

### Warranty Service

All Lifelong consumer products are covered against manufacturing defects from their date of purchase.

Name of the product: \_\_\_\_\_

Model: \_\_\_\_\_

Warranty coverage: 1 year

Please note: Purchase receipt is necessary for warranty verification.

Customer Care: [customercare@lifelongindia.com](mailto:customercare@lifelongindia.com)

### Customer Details

Name: \_\_\_\_\_

Address: \_\_\_\_\_

Home Number: \_\_\_\_\_

Office Number: \_\_\_\_\_

E-mail Address: \_\_\_\_\_

### Product Details

Model No.: \_\_\_\_\_

Serial No.: \_\_\_\_\_

Purchase Date: \_\_\_\_\_

Invoice Number: \_\_\_\_\_

Online Site: \_\_\_\_\_

Please log on to [www.lifelongindiaonline.com](http://www.lifelongindiaonline.com) and complete the online warranty form with your personal and product details within 14 days.

### Terms and Conditions:

- This warranty is void if:
  - The completely filled warranty card is not presented at the time of servicing the product.
  - The product is not operated according to the instructions given under the user manual.
  - Damages are caused by lightning, abnormal voltage, water or other liquid intrusion, fire, flood, accident, negligence or improper handling.
  - Product has been damaged due to installation, repairs, alterations or modifications by unauthorized service organizations or persons.
  - Product label specifying the model number, serial number and production code has been removed and altered.
  - Defects or parts requiring replacement due to ordinary wear and tear, corrosion, rust or stains, scratches, dents on the casing or paintwork of the product.
  - Claims damaged and/or missing parts (accessories) after 7 days from the original date of goods received.
  - Defects or faults in product which have been used for commercial/industrial purposes or which have been rented/leased or have been otherwise subject to non-household/non-domestic use.
- Repairs or replacements will be carried out by authorized service provider.
- During the limited warranty period, Lifelong or its authorized service provider will repair without charge the defective unit inclusive of labor and parts and restore the unit to its optimum working condition. All defective parts used for the warranty repair should be surrendered to Lifelong and/or its authorised servicer.
- All expenses incurred in collecting the unit (s) or part(s) thereof from authorised service provider as well as any other expenses and incidentals will be borne by the consumer.
- Lifelong obligations are limited to the repair and replacement of the defective product. Except as set forth above, there are no other express or implied warranty and all warranted, conditions or other terms implied by statute or common law (including any warranty of satisfactory quality, merchantability or fitness for a particular purpose) are excluded to the fullest extent permitted by the law.
- Lifelong total liability for damages relating to or arising out of the purchase or use of the product regardless of the type or cause of such damage of the form of characterization of the claim asserted (e.g. contract or tort) shall not exceed the original purchase price paid for the product.
- However in no event shall Lifelong and Lifelong authorized distributors be liable for any punitive, special incidental, indirect or consequential losses or damages whatsoever (including without limitation, damages for the loss of revenue, business, profits, goodwill, or contracts, business interruptions, loss of business information or any other pecuniary loss.) Whether or not Lifelong has been advised of the possibility of such damages. These limitations shall apply not with standing the failure of the essential purpose of any limited warranty. This limited warranty does not affect consumer's statutory rights under the law.
- No carrier, dealer or employee is thereof authorized to make modifications to this warranty and you should not rely on any such representation. Lifelong reserves the right to amend the terms and conditions if necessary.