

Manufactured for & Marketed by :  
**Lifelong Online Retail Private Limited**

For queries and complaints:  
please contact: [customercare@lifelongindia.com](mailto:customercare@lifelongindia.com)  
[www.lifelongindiaonline.com](http://www.lifelongindiaonline.com)

# lifelong

## FOOT AND CALF MASSAGER

### User Manual



**Built for Relaxation. Built for India.**

## INTRODUCTION:

Welcome to the Lifelong family.

Dear Customer,

Congratulations on purchasing Lifelong Foot & Calf Massager and a warm welcome to the Lifelong family.

At Lifelong, we take pride in developing quality products for your home, health and happiness. Through constant innovation in our products, we have been able to develop a wide spectrum of products across Home, Kitchen, Grooming and Wellness categories. We strive to enrich the lives of the modern consumer every day.

Keeping that in mind, we have offered all our products at a fair and honest price.

We hope that this Foot & Calf Massager serves you and your family's needs without any hassles.

Thanks again for becoming a part of the Lifelong family. Enjoy using the product!

Warm Regards,

## INSTRUCTIONS:

- Ensure that the device is switched off before you plug in.
  - Do not stand on, or exert excessive pressure on the massager, it can slow down or stop the machine
  - You may feel itchy while use. Do not worry as this is due to improved circulation
  - Do not use the massager while you're in the washroom
  - Do not use the massager on wet body, this may cause short circuit or electric shock
  - Those suffering from cancer, osteoporosis, tumors, high fever or pregnancy must consult a doctor before use.
  - Do not use the massager when it is in a bad condition
  - Never operate the device especially if there are children in the vicinity.
  - Do not use with animals or pets.
  - Do not cover the massager with the damp cloth when it is in an operating condition, it may result to a subsequent breakdown or fire.
  - Connect the device only to the type of voltage stated on the nameplate.
  - Do not sleep while using the massager
- As with every electrical device, this massage device should be used carefully and cautiously in order to avoid the risk of electric shock:
- Unplug the mains cable immediately after use.

Prevent any contact with water (except when cleaning with a slightly damp cloth). Water must never get inside the device. Never submerge the device in water. Never use the device in the bath, in the shower, in a swimming pool or over a basin filled with water. However, if water gets inside the casing, unplug the device immediately and contact customer service.

- The device and in particular the mains cable should be checked for any damage before each use.
- Never plug the device into the mains power supply if the cable is damaged.
- Never hold or carry the device by mains cable.
- If a fault should occur during use, switch off the device immediately and unplug it from the socket.
- The device is intended only for the purpose described in these instructions for use. The manufacturer can accept no liability for injury arising from incorrect or careless use.

## HOW TO USE:

The Lifelong Foot and Calf Massager can give you effective foot massage to your feet, ankles and calves all at the same time.

STEP 1: Ensure that the device is switched off before plug in.

STEP 2: Put your feet and calves in the black cloth area.

STEP 3: Switch on the device by sliding the switch to ON position

STEP 4: You can change the settings by clicking on this buttons.

STEP 5: It has 3 buttons to choose between Vibration & Heat, Manual and Auto.

STEP 6: The auto mode has 3 preset actions P1, P2 & P3

STEP 7: Switch to manual mode and adjust heat, vibration and speed on calves and feet.



## CLEANING:

- To clean the unplugged device after it has cooled down, wipe with a slightly dry cloth. However, ensure that water never gets inside the device.
- The foot cloth can be removed from the zipper and washed separately.

## SPECIFICATIONS:

Model No.	LLM99
Voltage	AC 220V-240V
Power	80W
Rated Working Time	15 Minutes

## Consumer Warranty Card

Dear Customer,

Thank you for choosing a Lifelong consumer product. All Lifelong Consumer products are designed and manufactured to the highest standards to deliver high quality performance, as well as easy installation and use. At Lifelong, we believe in providing not only service, but adding value to your purchase. The warranty has therefore been designed especially for you with your interest at heart.

### Warranty Service

All Lifelong consumer products are covered against manufacturing defects from their date of purchase.

Name of the product: \_\_\_\_\_

Model: \_\_\_\_\_

Warranty coverage: 1 year

Please note: Purchase receipt is necessary for warranty verification.

Customer Care: [customercare@lifelongindia.com](mailto:customercare@lifelongindia.com)

### Customer Details

Name: \_\_\_\_\_

Address: \_\_\_\_\_

Home Number: \_\_\_\_\_

Office Number: \_\_\_\_\_

E-mail Address: \_\_\_\_\_

### Product Details

Model No.: \_\_\_\_\_

Serial No.: \_\_\_\_\_

Purchase Date: \_\_\_\_\_

Invoice Number: \_\_\_\_\_

Online Site: \_\_\_\_\_

Please log on to [www.lifelongindiaonline.com](http://www.lifelongindiaonline.com) and complete the online warranty form with your personal and product details within 14 days.

### Terms and Conditions:

1. This warranty is void if:
  - a. The completely filled warranty card is not presented at the time of servicing the product.
  - b. The product is not operated according to the instructions given under the user manual.
  - c. Damages are caused by lightning, abnormal voltage, water or other liquid intrusion, fire, flood, accident, negligence or improper handling.
  - d. Product has been damaged due to installation, repairs, alterations or modifications by unauthorized service organizations or persons.
  - e. Product label specifying the model number, serial number and production code has been removed and altered.
  - f. Defects or parts requiring replacement due to ordinary wear and tear, corrosion, rust or stains, scratches, dents on the casing or paintwork of the product.
  - g. Claims damaged and/or missing parts (accessories) after 7 days from the original date of goods received.
  - h. Defects or faults in product which have been used for commercial/industrial purposes or which have been rented/leased or have been otherwise subject to non-household/non-domestic use.
2. Repairs or replacements will be carried out by authorized service provider.
3. During the limited warranty period, Lifelong or its authorized service provider will repair without charge the defective unit inclusive of labor and parts and restore the unit to its optimum working condition. All defective parts used for the warranty repair should be surrendered to Lifelong and/or its authorised servicer.
4. All expenses incurred in collecting the unit (s) or part(s) thereof from authorised service provider as well as any other expenses and incidentals will be borne by the consumer.
5. Lifelong obligations are limited to the repair and replacement of the defective product. Except as set forth above, there are no other express or implied warranty and all warranted, conditions or other terms implied by statute or common law (including any warranty of satisfactory quality, merchantability or fitness for a particular purpose) are excluded to the fullest extent permitted by the law.
6. Lifelong total liability for damages relating to or arising out of the purchase or use of the product regardless of the type or cause of such damage of the form of characterization of the claim asserted (e.g. contract or tort) shall not exceed the original purchase price paid for the product.
7. However in no event shall Lifelong and Lifelong authorized distributors be liable for any punitive, special incidental, indirect or consequential losses or damages whatsoever (including without limitation, damages for the loss of revenue, business, profits, goodwill, or contracts, business interruptions, loss of business information or any other pecuniary loss.) Whether or not Lifelong has been advised of the possibility of such damages. These limitations shall apply not withstanding the failure of the essential purpose of any limited warranty. This limited warranty does not affect consumer's statutory rights under the law.
8. No carrier, dealer or employee is thereof authorized to make modifications to this warranty and you should not rely on any such representation. Lifelong reserves the right to amend the terms and conditions if necessary.