


Lifelong

# Lifelong Mini Massager



Easy to  
hold & use



Unique pressure  
activated design



High performace  
vibrator

1 YEAR  
WARRANTY

## INTRODUCTION:

Welcome to the Lifelong family,

Dear Customer,

Congratulations on purchasing Lifelong Mini Massager and a warm welcome to the Lifelong family.

At Lifelong, we take pride in developing quality products for your home, health and happiness. Through constant innovation in our products, we have been able to develop a wide spectrum of products across Home, Kitchen, Grooming and Wellness categories. We strive to enrich the lives of the modern consumer every day. Keeping that in mind, we have offered all our products at a fair and honest price.

We hope that the Mini Massager serves you and your family's needs without any hassles.

Thanks again for becoming a part of the Lifelong family. Enjoy using the product!

Warm Regards,

Team Lifelong

## PRODUCT PRINCIPAL:

- Enjoy effortless relaxing massage anywhere, any time with this vibrating mini massager.
- Its ergonomic design and its 3 robotic massage contacts allow a proper massage all along your body. Also, its vibrating movement benefits in blood circulation and helps in muscle relaxation. When combined with a balanced diet and exercise, the massage helps to decrease fat deposits.
- Recommended massage session may last from 15 minutes to 1 hour.



## PRODUCT STRUCTURE:

This machine is well designed based on contour lines of the body with three massage fingers and a high performance motor, which generates suitable and wanted massage effects controlled by a computer chip.

**Features:** The instrument is designed in waterproof structure. The whole machine is waterproof.

## KNOW YOUR PRODUCT:



1. Main Unit



2. Batteries needed  
(AAAx3)



3. User Manual


## METHOD OF OPERATION:

### BATTERY INSTALLATION

1. Open the cover of the main machine according to the indicated direction and put 3 AAA batteries inside. (Battery not included)

2. Close the battery cover

- Please confirm the anode and cathode of batteries.
- Please refer to above chart as for anode and cathode of batteries.



### OPERATION:

1. Put the 3 massage fingers of the massager and press the top of the main machine with proper force, then it start working.

2. Remove this machine from your body without any pressure on its top, then it ceases working.

## TRANSPORATION AND STORAGE:

### Transportation:

Ambient temperature: 5 °C-40

Relative humidity: not more than 80%, ventilated room and without corrosive gas.

Atmosphere pressure: 860hpa-1060hpa

This machine is equipped with simple shockproof facility, available for airplane, railway, road and ship transport. It should avoid splash, inversion and collision.

## TROUBLESHOOTING:

Fault	Reason	Solution
Not Working	Low batteries? Wrong direction of the anode and cathode of battery?	Change the batteries correct direction of the anode and cathode of batteries.
Weak Massage	Low batteries?	Change batteries.

## MAINTENANCE:

- Moisten a soft cloth with medical alcohol to wipe the machine body to remove the dirt.
- Please keep this machine on safe place after use. No heavy weight is allowed on the machine.
- Please ask your professional technician for repair if any fault occurs. Do not disassemble the instrument without permission.

## PRODUCT SPECIFICATIONS:

- Dimension: 105.6mm x 6.4mm x 89.8mm
- Battery: AAA dry battery (Battery not included)
- Net Weight: 204.5
- Materials: ABS resin gum
- Power consumption: 1.3 VA
- Ingression Protection Rating: IP6 7

## Consumer Warranty Card

Dear Customer,

Thank you for choosing a Lifelong consumer product. All Lifelong Consumer products are designed and manufactured to the highest standards to deliver high quality performance, as well as easy installation and use. At Lifelong, we believe in providing not only service, but adding value to your purchase. The warranty has therefore been designed especially for you with your interest at heart.

### Warranty Service

All Lifelong consumer products are covered against manufacturing defects from them date of purchase.

Name of the product: \_\_\_\_\_

Model: \_\_\_\_\_

Warranty coverage: 1 year

Please note: Purchase receipt is necessary for warranty verification.

Customer Care: [customercare@lifelongindia.com](mailto:customercare@lifelongindia.com)

### Customer Details

Name: \_\_\_\_\_

Address: \_\_\_\_\_

Home Number: \_\_\_\_\_

Office Number: \_\_\_\_\_

E-mail Address: \_\_\_\_\_

### Product Details

Model No.: \_\_\_\_\_

Serial No.: \_\_\_\_\_

Purchase Date: \_\_\_\_\_

Invoice Number: \_\_\_\_\_

Online Site: \_\_\_\_\_

Please log on to [www.lifelongindiaonline.com](http://www.lifelongindiaonline.com) and complete the online warranty form with your personal and product details within 14 days.

### Terms and Conditions:

1. This warranty is void if:
  - a. The completely filled warranty card is not presented at the time of servicing the product.
  - b. The product is not operated according to the instructions given under the user manual.
  - c. Damages are caused by lightning, abnormal voltage, water or other liquid intrusion, fire, flood, accident, negligence or improper handling.
  - d. Product has been damaged due to installation, repairs, alterations or modifications by unauthorized service organizations or persons.
  - e. Product label specifying the model number, serial number and production code has been removed and altered.
  - f. Defects or parts requiring replacement due to ordinary wear and tear, corrosion, rust or stains, scratches, dents on the casing or paintwork of the product.
  - g. Claims damaged and/or missing parts (accessories) after 7 days from the original date of goods received.
  - h. Defects or faults in product which have been used for commercial/industrial purposes or which have been rented/leased or have been otherwise subject to non-household/non-domestic use.
2. Repairs or replacements will be carried out by authorized service provider.
3. During the limited warranty period, Lifelong or its authorized service provider will repair without charge the defective unit inclusive of labor and parts and restore the unit to its optimum working condition. All defective parts used for the warranty repair should be surrendered to Lifelong and/or its authorised servicer.
4. All expenses incurred in collecting the unit (s) or part(s) thereof from authorised service provider as will as any other expenses and incidentals will be borne by the consumer.
5. Lifelong obligations are limited to the repair and replacement of the defective product. Except as set forth above, there are no other express or implied warranty and all warranted, conditions or other terms implied by statute or common law (including any warranty of satisfactory quality, merchantability or fitness for a particular purpose) are excluded to the fullest extent permitted by the law.
6. Lifelong total liability for damages relating to or arising out of the purchase or use of the product regardless of the type or cause of such damage of the form of characterization of the claim asserted (e.g. contract or tort) shall not exceed the original purchase price paid for the product.
7. However in no event shall Lifelong and Lifelong authorized distributors be liable for any punitive, special incidental, indirect or consequential losses or damages whatsoever (including without limitation, damages for the loss of revenue, business, profits, goodwill, or contracts, business interruptions, loss of business information or any other pecuniary loss.) Whether or not Lifelong has been advised of the possibility of such damages. These limitations shall apply not with standing the failure of the essential purpose of any limited warranty. This limited warranty does not affect consumer's statutory rights under the law.
8. No carrier, dealer or employee is thereof authorized to make modifications to this warranty and you should not rely on any such representation. Lifelong reserves the right to amend the terms and conditions if necessary.

Manufactured for & Marketed by :  
**Lifelong Online Retail Private Limited**

For queries and complaints:  
please contact: [customercare@lifelongindia.com](mailto:customercare@lifelongindia.com)  
[www.lifelongindiaonline.com](http://www.lifelongindiaonline.com)