FACIAL STEAMER & VAPORIZER

User Manual

 ${\it Manufactured for \& Marketed by:}$

Lifelong Online Retail Private Limited

For queries and complaints: please contact: customercare@lifelongindia.com www.lifelongindiaonline.com



Built for Relaxation. Built for India.

INTRODUCTION:

Dear Customer,

Congratulations! We welcome you to the Lifelong family, as you welcome our Lifelong Facial Steamer & Vaporizer to be a part of your home.

At Lifelong, our products are carefully built for the needs of India, which is why we have crafted a diverse range of innovative products, across a spectrum of categories - Home, Kitchen, Grooming, Fitness, Lifestyle and Smart Home, to better help you.


We then deliver them at the best prices, enriching your everyday living.

We are sure you can't wait to start using your new Lifelong Facial Steamer & Vaporizer.

We know we can't!

Looking forward to building our bond for life, cheers!

Thank you, Team Lifelong



LLS36

1

TABLE OF CONTENTS:

- 1. About the product
- 2. Precautions Before Use
- 3. Operating Instructions
- 4. Device Maintenance
- 5. Terms & Conditions
- 6. Warranty Card

ABOUT THE PRODUCT:

Lifelong Facial Steamer & Vaporizer is designed for steam inhalation to provide relief from cold and cough, respiratory problems, sinusitis, and nasal congestion, Apart from this, it can also be used for steaming face which helps in removing impurities, excess oil from the face, penetrating into your skin to remove blackheads and unclog pores.

PRECAUTIONS BEFORE USE:

- If you have heart expansion disease, serious rosacea, eczema, or psoriasis, please do not use this product.
- · Do not add salt into the water
- Pregnant women and patients with serious diseases, such as in the use of the products need to add herbs and essential oil before, please consult a physician.
- Once appear the skin inflammation or discomfort, immediately stop using essential oils and herbs.
- Often clean and detergent can ensure the best effect and extend the product life
- · Always unplug the steamer immediately after use.
- Do not place in or drop into water or other liquids. Never use while bathing.
- Do not reach for an appliance that has fallen into water or other liquids. Firstly, unplug it immediately.

OPERATING INSTRUCTIONS:

• Fill the Water Tank:

The steamer comes with a water tank. Fill the tank with lukewarm distilled water instead of regular tap water.

• Turn on the Steamer:

Allow the water to heat up so that steam starts coming out from the facial steamer.

Steam Your Face:

Allow the steam to fall over your face from a certain distance. You can control the steam delivery as per your liking by adjusting the control knob.

• Steam Timing:

Do this for about 3 to 10 minutes, depending on your skin's type, needs and reaction.

 Do not use it for more than 5 minutes if you have sensitive skin or have any skin problems.

DEVICE MAINTENANCE:

- After each use, it should be cleaned properly to avoid oil and other dirt deposition.
- Always turn off and unplug it from the electrical outlet before you start cleaning.
- Let it cool down, then empty excess water, and wipe dry interior of the facial mask and rim around the opening.

Terms and Conditions:

1.This warranty is void if:

- a. The completely filled warranty card is not presented at the time of servicing the product.
- The product is not operated according to the instructions given under the user manual.
- C. Damages are caused by lightning, abnormal voltage, water or other liquid intrusion, fire, flood, accident, negligence or improper handling.
- d. Product has been damaged due to installation, repairs, alterations or modifications by unauthorized service organizations or persons.
- e. Product label specifying the model number, serial number and production code has been removed and altered.
- Defects or parts requiring replacement due to ordinary wear and tear, corrosion, rust or stains, scratches, dents on the casing or paintwork of the product.
- g. Claims damaged and/or missing parts (accessories) after 7 days from the original date of goods received.
- Defects or faults in product which have been used for commercial/industrial purposes or which have been rented/leased or have been otherwise subject to non-household/non-domestic use.
- Repairs or replacements will be carried out by authorized service provider.
- During the limited warranty period, Lifelong or its authorized service provider will repair without charge the defective unit inclusive of labor and parts and restore the unit to its optimum working condition. All defective parts used for the warranty repair should be surrendered to Lifelona and/or its authorised servicer.
- All expenses incurred in collecting the unit (s) or part(s) thereof from authorised service provider as will as any other expenses and incidentals will be

- 5. Lifelong obligations are limited to the repair and replacement of the defective product. Except as set forth above, there are no other express or implied warranty and all warrantied, conditions or other terms implied by statute or common law (including any warranty of satisfactory quality, merchantability or fitness for a particular purpose) are excluded to the fullest extent permitted by the lass.
- 6. Lifelong total liability for damages relating to or arising out of the purchase or use of the product regardless of the type or cause of such damage of the form of characterization of the claim asserted (e.g. contract or tort) shall not exceed the original purchase price paid for the product.
- 7. However in no event shall Lifelong and Lifelong authorized distributors be liable for any punitive, special incidental, indirect or consequential losses or damages whatsoever (including without limitation, damages for the loss of revenue, business, profits, goodwill, or contracts, business interruptions, loss of business information or any other pecuniary loss). Whether or not Lifelong has been advised of the possibility of such damages. These limitations shall apply not with standing the failure of the essential purpose of any limited warranty. This limited warranty does not affect consumer's statutory rights under the law.
- No carrier, dealer or employee is thereof authorized to make modifications to this warranty and you should not reply on any such representation. Lifelong reserves the right to amend the terms and conditions if necessary.

6

Consumer Warranty Card

Dear Customer,

Thank you for choosing a Lifelong consumer product. All Lifelong Consumer products are designed and manufactured to the highest standards to deliver high quality performance, as well as easy installation and use. At Lifelong, we believe in providing not only service, but adding value to your purchase. The warranty has therefore been designed especially for you with your interest at heart.

date of purchase. Name of the produc	ner products are covered against manufacturing defects from them
Warranty coverage	
	se receipt is necessary for warranty verification. tomercare@lifelongindia.com
Address: Home Number: Office Number:	
Serial No.: ————————————————————————————————————	
Purchase Date: — Invoice Number: — Online Site: —	

Please log on to www.lifelongindiaonline.com and complete the online warranty form with your personal and product details within 14 days.